3M ESPE

Ketac™ Molar Aplicap™

Glass Ionomer Restorative Material

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ENGLISH

Product Description

Ketac™ Molar Aplicap™ is a radiopaque bulk fill glass ionomer cement in capsules. On account of its chemical adhesion to the dental enamel and dentin, a gentle preparation and an especially tight seal at the margins is possible. The product is available in various shades; it can be applied without lining and releases fluoride ions. The dispensable quantity of a capsule is at least 0.1 ml.

For these Instructions for Use should be kept for the duration of product use. For details on all mentioned products, please refer to the respective Instructions for Use. Ketac™ Glaze, Scotchbond™ Universal Adhesive, or Single Bond Universal Adhesive are not available in all countries.

Indications

- Linings for single-surface and multiple-surface composite fillings
- Core build-up prior to crown placement
- Primary tooth fillings
- Single-surface fillings in non-occlusal load-bearing areas
- Cervical fillings. If aesthetics is not the prime consideration
- Uniplanar and multi-planar temporary fillings

Precautionary Measures

3M MSDSs can be obtained from www.mmm.com, or contact your local subsidiary.

Cavity Preparation

- Remove carious tooth substance only; undercuts are not required.
- Do not prepare thin beveling margins. A wall thickness of 0.5 mm has to be provided for the material.

- In order to prevent pulpal irritation, cover areas in close proximity to the pulp with local application of a calcium hydroxide material if necessary.

Conditioning

- For optimal chemical bonding to the teeth, the smear layer produced in the preparatory steps must be carefully removed by applying Ketac™ Conditioner to the prepared surfaces and allowing the substance to react for 10 sec.
- Then rinse with copious amounts of water.
- Blow the cavity dry in only 2-3 short intervals with water- and oil-free air or dab dry with cotton pellets. Do not overdry! The cavity should have a matt shiny appearance. Excessive drying can result in post-operative sensitivity after filling.
- Avoid further contamination.

Times

The following times apply for an ambient temperature of 23° C/73° F:

	min:sec
Activation	0:02
Mixing in RotoMix™	0:10
Mixing in high-frequency mixer, e.g., CapMix™	0:15
Working after start of mixing	2:00
Setting after start of mixing	4:30

Setting is accelerated at room temperatures over 23° C/73° F. Exceeding the given working time will cause the loss of adhesion to the dental enamel and the dentin.

Capsule Activation

- Place the Aplicap™ Activator on a sturdy surface and insert the Aplicap capsule into the activator. In doing so, keep the applicator nozzle shut.
 - Both the Aplicap Activator and the Applier are color-coded in orange.
- Using the ball of your hands, depress the activator lever firmly and completely to the stop and hold it down for 2 to 4 seconds.
 - Firmly and completely depressing the lever to the stop and holding it in this position is the only way to ensure that all of the liquid is pressed into the powder. If inadequate pressure is applied or if the lever is not depressed completely, it is possible that too little liquid will be mixed with the powder, increasing the viscosity and changing the product characteristics.

Mixing

- Mix the capsule at approx. 4,300 rpm in a high-frequency mixing device such as CapMix™ or the RotoMix rotating mixer. In doing so, keep the application nozzle shut.

Application

- Insert the capsule into the Aplicap™ Applier and open the application nozzle.
- Apply the glass ionomer cement in the cavity and shape the filling.
 - Ensure that there is no contamination from water, saliva, etc., for the entire duration of application and shaping; the working environment must be kept dry.

Attention

- By all means prevent delay between activation, mixing, and application because the beginning of setting of the material in the capsule will render difficult or prevent the extrusion.
- Ketac Molar Aplicap clings to metal instruments and for this reason should be washed off with cold water before setting.

Filling Protection

Glass ionomers go through two phases of setting and are susceptible to water during the first setting period. Subsequent dehydration may be possible. Hence extra care should be taken to avoid moisture exposure or trying out, especially at the end of the working time. In the cases where this is not possible, Ketac Glaze or Scotchbond Universal Adhesive or Single Bond Universal Adhesive is used as filling protection.

- Immediately after shaping, apply Ketac Glaze or Scotchbond Universal Adhesive or Single Bond Universal Adhesive to all free surfaces of the glass ionomer cement and light-cure according to the respective Instructions for Use.
- The inhibition layer remaining after the setting of Ketac Glaze or Scotchbond Universal Adhesive or Single Bond Universal Adhesive may obstruct the setting of the addition-type silicone impression materials. For this reason, remove this layer or use a polyether impression material if needed.

Finishing

- Further finishing with Arkansas stones, fine-grain diamonds, abrasive discs with successively decreasing grain size, or silicone polishers should be started 4:30 minutes after start of mixing at the earliest.
- If necessary, apply Ketac Glaze or Scotchbond Universal Adhesive or Single Bond Universal Adhesive again for surface sealing and light-cure according to the respective Instructions for Use.

Storage and Shelflife

Store the capsules in the blister packs. Use within 1 month maximum once the blister pack has been opened! Store the product at 15–25° C/59–77° F. Do not use after the expiry date.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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