3M ESPE

Photac™ Fil Quick Aplicap™

Lichthärtendes Glasionomer Dental Restorative

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3M ESPE, Aplicap, CapMix, Elipar, Ketac, Photac and RotoMix are trademarks of 3M or 3M ESPE AG.
© 3M 2011. All rights reserved.

ENGLISH

Product description

Photac™ Fil Quick Aplicap™, manufactured by 3M ESPE, is a dual-plus curing glass ionomer cement in a capsule. Since it bonds chemically to enamel and dentin, it permits gentle preparation. It is available in eight shades. The material can be applied without cavity lining. It releases fluoride ions and is radiopaque. The applicable quantity of a capsule amounts to at least 0.1 ml. For details on all products mentioned below please refer to the corresponding instructions for Use.

☞ Instructions for Use should not be discarded for the duration of the product use.

Indications

- Fillings in Black's class III and V
- Small fillings in Black's class I
- Wedge-shaped defects
- Fillings in deciduous teeth
- Fissure sealing
- Temporary restorations in Black's class I and II

Precautionary Measures

For Patients

Photac Fil Quick Aplicap liquid contains substances that may cause allergic reactions by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies.

If prolonged contact with oral tissue occurs, flush with large amounts of water. If an allergic reaction occurs, remove the product, and discontinue future use; seek medical attention as needed.

For Dental Personnel

Photac Fil Quick Aplicap liquid contains substances that may cause allergic reactions by skin contact in certain individuals. To reduce the risk of allergic response, avoid exposure to these materials; in particular, avoid exposure to uncured materials. If skin contact occurs, wash skin with soap and water. The use of protective gloves and a no-touch technique is recommended.

Acrylates may penetrate commonly used gloves. If the liquid contacts the glove, remove and discard the gloves. Wash the hands immediately with soap and water and re-glove.

If an allergic reaction occurs, seek medical attention as needed.

3M ESPE MSDSs can be obtained at www.mmm.com or from your local subsidiary.

- Remove caries dentin only. Undercuts are not required.
- Use water spray to remove debris.
- Blow the cavity dry in only 2-3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry! **The cavity should be just dry but the surface has a matt shiny appearance.** Excessive drying can result in post-operative sensitivity after filling.
- Avoid any contamination by saliva.

Pulp Protection

► To prevent pulpal irritation, the areas around the pulp should be covered by dabbing them with a hardsetting calcium hydroxide preparation.

Conditioning

► No further preparatory procedures are required. In case of an unfavorable cavity geometry, the use of Ketac™ Conditioner, manufactured by 3M ESPE, can further improve the adhesive stability obtainable.

Activation

- Place the Aplicap™ Activator, manufactured for 3M ESPE, on a sturdy surface and insert the Aplicap capsule in the activator.
- Using the ball of your hands, depress the activator lever firmly and completely to the stop and hold it down for 2 to 4 seconds:
 - Firmly and completely depressing the lever to the stop and holding it in this position is the only way to ensure that all of the liquid is pressed into the powder. If inadequate pressure is applied or if the lever is not depressed completely, it is possible that too little liquid will be mixed with the powder, increasing the viscosity and changing the product characteristics.

Mixing
► Mix the capsules at approx. 4.300 rpm, in a high frequency mixing device such as CapMix™ (for 10 sec) or the RotoMix™ rotating mixer (for 8 sec), and apply with the Aplicap™ Applier, all products manufactured by or for 3M ESPE. Both the Aplicap Activator and the Aplicap Applier are color-coded in red.

Warning!
It is imperative to avoid waiting times between activation, mixing, and application as the onset of the curing prevents or renders the extrusion of the material difficult.

- Use Aplicap™ applier to **apply layers of up to 2 mm maximum**. Avoid all types of contamination during application.
- Use halogen or LED light (such as Elipar™ FreeLight, manufactured by 3M ESPE) for curing.
 - Place the light source as closely as possible to the filling.
 - The maximum layer thickness (without intermediate curing) may not exceed 2 mm.
 - In case of cavities of more than 2 mm in depth, several layers have to be applied.
 - The self-curing mechanism ensures that the material also cures in areas sheltered from the light. In order to attain its optimum physical properties and to ensure immediate preparation, the material has to be light-cured carefully.
 - It is possible to add layers on material already cured if the surface is clean and dry.
- The material adheres to metal instruments and thus has to be washed off in cold water before it sets.

Times

The following times apply at a room temperature of 23° C/73° F:

	min.:sec.
Activation	0:02
Mixing in RotoMix	0:10
Mixing in a high-frequency mixer, such as CapMix	0:15
Manipulation from the beginning of the mixture	2:00
Light curing	0:20 (DBO 0:40)

At room temperatures of over 23° C/73° F the curing is accelerated. Shorter mixing times prolong the manipulation and curing time, whereas longer mixing times shorten the manipulation and curing time. The mixing time in high-frequency mixers can be selected at between 8 and 15 seconds, and the mixing time in RotoMix can be set to between 7 and 12 seconds. If the manipulation time is exceeded, bonding to the enamel and dentin is lost.

Finishing

- Use a sharp instrument to remove excess material.
- Use rotary instruments (such as polishing discs or silicone polishers) to finish under cold water.

Note

- Humidity in the cavity prevents bonding significantly.
- Photac Fil Quick Aplicap polymerises in daylight or artificial lighting. Therefore, the material should not be dispensed from the capsule until immediately prior to manipulation. **The manipulation time is reduced drastically in case of application under the operating light!**

Storage and shelf-life

Store the product at 15–25° C/59–77° F. Do not use after the expiration date.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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