Variolink® Esthetic is a colour-stable, adhesive luting system for the permanent cementation of ceramic and composite resin restorations. Variolink Esthetic is offered in a purely light-curing (Variolink Esthetic LC) and a dual-curing version (Variolink Esthetic DC). The special filler composition gives Variolink Esthetic a very high radiopacity.

Variolink Esthetic is available in the following five shade gradations:

Shade		Effect
Light +	Translucency approx. 5% (opaque)	Lighter / more whitish
Light	Translucency approx. 10%	
Neutral	Translucency approx. 17%	None
Warm	Translucency approx. 12%	Darker / more yellowish
Warm +	Translucency approx. 8.5% (opaque)	

The gradation of Variolink Esthetic shades is based on the effect that a luting composite exerts on the brightness value of the final restoration. Variolink Esthetic Neutral does not affect the brightness value. At the same time, it demonstrates the highest translucency and is thus neutral in shade. Light and Light + make the restorations lighter, while Warm and Warm + create a darker overall shade impression. In order to evaluate the overall effect of the restoration in conjunction with the various Variolink Esthetic shades prior to permanent cementation, we recommend using the Variolink Esthetic try-in pastes.

The working and curing times depend on the ambient temperature. The following times apply as soon as Variolink Esthetic DC has been extruded from the automix syringe:

When used purely self-curing	Room temperature 23°C ± 1°C/ 73°F ± 2°F	Intraoral 37°C ± 1°C/ 99°F ± 2°F
Working time	approx. 5 min	approx. 2 min
Setting time (incl. working time)	approx. 8 min	approx. 4 min

Variolink Esthetic DC is always dispensed from the automix syringe in the optimum ratio.

Variolink® Esthetic DC

EN Instructions for Use	- Dual-curing resin-based dental luting material	DA Brugsanvisning	- dualhærdende permanent plastrcement	SL Za navodila za uporabo običite	SQ Udhëzimet e përdorimit do t'i gjeni në
DE Gebrauchsinformation	- Dualhärtendes zahnärztliches Befestigungskomposit	FI Käyttöohjeet	- Kaksioskovaiteinen hartspohjainen hammaskilnitysmateriaali	HR Upute za upotrebu možete pronaći na adresi:	RO Pentru instrucțiuni de utilizare, vă rugăm să utilizați link-ul de mai jos
FR Mode d'emploi	- Composite de collage à prise duale à usage dentaire	NO Bruksanvisning	- Dobbeltstående, harpiksbasert dentalt festemateriale	CS Návod k použití naleznete na internetové stránce	UK Для ознайомлення з інструкціями із застосування, будь ласка, відвідайте сайт
IT Istruzioni d'uso	- Cemento composito dentale ad indurimento duale	NL Productinformatie	- Tandheelkundig, dual uithardend bevestigingsmateriaal op basis van kunsthar	HU Kérjük, használati utasításokért keresse fel a	ET Kasutusjuhend asub aadressil
ES Instrucciones de uso	- Material dental de cementado en base a resina de polimerización dual	EL Οδηγίες Χρήσεως	- Οδοντιατρική ρητίνωδης κολλητική πολυμερισμού	SR Упутство за употребу можете наћи на интернет адреси	LV Lietošanas instrukcijas skatiet vietnē
PT Instruções de Uso	- Material para cimentação à base de resina de cura-dual	TR Kullanma Talimatı	- Dual sertleşen rezin esaslı dental	MK Упатството за употреба се наоѓа на веб-страницата	LT Naudojimo instrukcijas rasite apsilankę

Composition

The monomer matrix of Variolink Esthetic is composed of urethane dimethacrylate and further methacrylate monomers. The inorganic fillers are ytterbium trifluoride and spheroid mixed oxide. Initiators, stabilizers and pigments are additional ingredients. The particle size is 0.04–0.2 µm. The mean particle size is 0.1 µm. The total volume of inorganic fillers is approx. 38%.

Indication

- Adhesive luting of glass-ceramic, lithium disilicate glass-ceramic and composite restorations (inlays, onlays, partial crowns, crowns, bridges).
- Restorations made of opaque ceramics, e.g. oxide ceramics, can only be permanently cemented if an adhesive is additionally used that is separately light-cured.
- The following overview provides recommendations for the selection of the suitable Variolink Esthetic shade:

		Indications	
Shade	Effect	Inlays / Onlays	Crowns / Bridges
Light +	considerably lightens	---*	✓
Light	slightly lightens	✓	✓
Neutral	no shade effect	✓	✓
Warm	slightly darkens/ more yellowish	✓	✓
Warm +	considerably darkens/ more yellowish	---*	✓

* The use of the shades Light+ and Warm+ may lead to visible restoration margins when cementing inlays/onlays.

Contraindication

The application of Variolink Esthetic is contraindicated

- if a dry working field cannot be established or the stipulated working procedures cannot be applied;
- if a patient is known to be allergic to any of the ingredients of Variolink Esthetic.

Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to individual components have been reported.

Interactions

Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the application of products containing these components, e.g. mouth rinses and temporary cements, must be avoided.

Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Therefore, do not disinfect the preparation and the syringe using oxidative agents. Disinfection may be carried out, for instance, by wiping with medical alcohol.

Alkaline jet mediums applied on dentin (e.g. Airflow) may compromise the effect of self-etching adhesives.

Application

For more detailed information also refer to the Instructions for Use of the products used in conjunction with Variolink Esthetic.

- 1 Removal of temporary restoration and cleaning of the cavity**

Remove possible residue of the temporary luting composite from the cavity or the preparation with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxyl[®] fluoride-free prophyl paste). Rinse with water spray.

Subsequently, dry with oil- and moisture-free air. Avoid overdrying.

Note:

Cleaning with alcohol can lead to dehydration of dentin.
- 2 Try-in of the restoration and isolation**

Seat the restoration using the desired Variolink Esthetic try-in paste and check the shade, fit and occlusion of the restoration.

Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture.

If necessary, make adjustments with fine diamonds at medium speed and with slight pressure and adequate water cooling. Polish ground surfaces.

When an adhesive luting protocol with composites is used, safe isolation of the operating field – preferably with a rubber dam, e.g. OptraDam[®] Plus, or alternatively with cotton rolls and a saliva ejector – is required. If the dental hard tissue is contaminated with saliva or blood during try-in, it needs to be cleaned again as described in section 1.
- 3 Pre-treatment of the restoration**
 - 3.1** If a restoration which has been pre-treated in the dental laboratory is contaminated with saliva or blood during try-in, the contaminated bonding surface of the restoration needs to be cleaned as follows:
 - Following try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
 - Shake Ivoclean before use and completely coat the bonding surface of the restoration with the product using a microbrush or bush.
 - Leave Ivoclean to react for 20 s, then thoroughly rinse with water spray and dry with oil-free air.
 - Next, prime the bonding surface of the restoration with a suitable bonding agent (e.g. Monobond[®] Plus).
Make sure to observe the Instructions for Use of the bonding agent used.
 - 3.2** Restorations that have not been pre-treated by the dental laboratory should be rinsed with water spray and dried after the try-in procedure.
 - 3.2.1** Subsequently, condition the bonding surface of the restoration as follows:
 - a) Glass-ceramic restorations (e.g. IPS Empress[®])**
 - Etch with 5% hydrofluoric acid (e.g. IPS[®] Ceramic Etching Gel) for 60 s or according to the instructions of the manufacturer of the restorative material.
 - Thoroughly rinse with water spray and dry with oil-free air.

- b) Lithium disilicate glass-ceramic restorations (e.g. IPS e.max® Press, IPS e.max CAD)
 - Etch with 5% hydrofluoric acid (e.g. IPS Ceramic Etching Gel) for 20 s or according to the instructions of the manufacturer of the restorative material.
 - Thoroughly rinse with water spray and dry with oil-free air.
- c) Zirconium oxide (e.g. IPS e.max ZirCAD, Zenostar®) or aluminium oxide ceramics
 - Sandblast the internal surface of the restoration (sandblasting parameters according to the instructions of the manufacturer of the restorative material).
 - If necessary, clean the restoration in an ultrasonic unit for about one minute.
 - Thoroughly rinse with water spray and dry with oil-free air.
 - **IMPORTANT!** In order to create a strong bond, do not clean the zirconium oxide surfaces with phosphoric acid.
- d) Composite or fibre-reinforced composite restorations
 - Blast the internal surface of the restoration (according to the instructions of the manufacturer of the restorative material).
 - If necessary, clean the restoration in an ultrasonic unit for about one minute.
 - Thoroughly rinse with water spray and dry with oil-free air.

3.2.2 Apply Monobond Plus with a brush or a microbrush to the pre-treated surfaces, let it react for 60 s and then disperse it with a strong stream of air.

4 Pre-treatment of the preparation and application of the adhesive (e.g. Adhese® Universal)

Please observe the Instructions for Use of the adhesive used.

4.1 Conditioning with phosphoric acid gel (optional)

The bond to enamel can be improved by selectively etching the enamel or by applying the "etch & rinse" technique. Unprepared enamel surfaces must be conditioned with phosphoric acid (e.g. Total Etch). Please observe the Instructions for Use for the phosphoric acid gel.

a) Selective enamel etching

Apply phosphoric acid gel (e.g. Total Etch) onto the enamel and allow it to react for 15–30 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

b) Etch & rinse technique

Apply phosphoric acid gel (e.g. Total Etch) onto the prepared enamel first, and then to the dentin. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

4.2 Application of Adhese Universal

- Starting with the enamel, thoroughly coat the tooth surfaces to be treated with Adhese Universal.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.

- Disperse Adhesive Universal with oil- and moisture-free compressed air until a glossy, immobile film layer results.
- Light-cure Adhese Universal for 10 s using a light intensity of $\geq 500 \text{ mW/cm}^2$ (e.g. Bluephase® Style).

5 Seating of the restoration

- For each application place a new mixing tip on the syringe. Dispense Variolink Esthetic DC from the automix syringe and apply the desired quantity directly to the restoration.
As the luting material will cure in the used mixing tip, it serves as a seal for the remaining contents of the syringe until needed again (replace with a new tip before the next use).
Note: Variolink Esthetic DC should be processed quickly after extrusion from the automix syringe and the restoration seated into place. Apply Variolink Esthetic DC directly to the internal surface of the restoration. The direct application of Variolink Esthetic DC to the preparation or cavity pretreated with adhesive is contraindicated, as this would significantly accelerate the curing process and could thus detrimentally affect the accuracy of fit of the restoration.
- Seat the restoration and **fix/hold it in place during excess removal**.
- Remove all excess luting material.
 - a) Wiping technique
Remove excess material immediately with a brush, dental floss or scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics).
 - b) Quarter technique (max. 6 bridge abutments) - Light-curing of excess and subsequent removal
Light-cure excess material with the polymerization light (e.g. Bluephase Style) **for 2 s** per quarter surface (mesio- oral, disto-oral, mesio-buccal, disto-buccal) at a distance of max. 10 mm.
Thereafter, excess cement is easy to remove with a scaler. Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal areas, gingival margins, pontics).

6 Polymerization

- As with all composite systems, Variolink Esthetic is subject to oxygen inhibition. This means that the surface layer (approx. 50 μm) does not polymerize during curing, as it comes in contact with atmospheric oxygen. To prevent this, cover the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after excess removal.
- Polymerize Variolink Esthetic in segments, starting with the proximal margins:

Light intensity	Exposure time per mm ceramic and segment	
500 - 1,000 mW/cm^2	20 s	
$\geq 1,000 \text{ mW/cm}^2$	10 s	e.g. Bluephase Style

- Rinse off glycerine gel/airblock (e.g. Liquid Strip).

7 Finishing of the completed restoration

- Finish the cement lines with preparation diamonds if necessary.
- Polish restoration margins/cement lines with suitable polishers (e.g. Astropol®) as well as finishing and polishing strips.
- Check occlusion and functions and make corrections if necessary.

Warning

- Prevent any contact of uncured Variolink Esthetic with the skin / mucous membrane and eyes.
- Uncured Variolink Esthetic may have a slight irritating effect and cause sensitization to methacrylates.
- Customary medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Store Variolink Esthetic at 2–28 °C/36–82 °F.
- Do not use Variolink Esthetic after the expiry date.
- Do not disinfect syringes with oxidizing disinfectants.
- In order to ensure a tight seal of Variolink Esthetic DC automix syringe, the mixing tip is left on the syringe after use.
- Expiry date: see note on syringes and packages.

Keep out of the reach of children!

For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

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www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai

Trading Co., Ltd., 2/F Building 1, 881 Wuding Road, Jing An District

200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 3399, Fax +57 1 633 1663

www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118, F-74410 Saint-Jorioz, France

Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, D-73479 Ellwangen, Jagst, Germany

Tel. +49 7961 889 0, Fax +49 7961 6326

www.ivoclarvivadent.de

Wieland Dental + Technik

GmbH & Co. KG, Schwenninger Strasse 13, D-75179 Pforzheim, Germany

Tel. +49 7231 3705 0, Fax +49 7231 3579 59

www.wieland-dental.com

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate

Veera Desai Road, A, ndheri (West), Mumbai, 400 053, India

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www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale, PO Box 303011 North Harbour

Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/, office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204

P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 11 293 8345, Fax +966 11 293 8344

www.ivoclarvivadent.com

Ivoclar Vivadent Pte. Ltd.

171 Chin Swee Road, #02-01 San Centre, Singapore 169877

Tel. +65 6535 6775, Fax +65 6535 4991

www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

C/ Ribera del Loira nº 46, 5ª planta, 28042 Madrid, Spain

Tel. + 34 913 757 820, Fax + 34 913 757 838

www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14, S-169 56 Solna, Sweden

Tel. +46 8 514 939 30, Fax +46 8 514 939 40

www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi, Sakayik Sokak, Nisantas' Plaza No:38/2

Kat:5 Daire:24, 34021 Sisli – Istanbul, Turkey

Tel. +90 212 343 0802, Fax +90 212 343 0842

www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building, Feldspar Close, Warrens Business Park,
Enderby, Leicester LE19 4SE, United Kingdom

Tel. +44 116 284 7880, Fax +44 116 284 7881

www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

www.ivoclarvivadent.us