G-aenial® Universal Injectable
UNIVERSAL LIGHT-CURED RADIOPAQUE INJECTABLE COMPOSITE
For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Direct restorative for Class I, II, III, IV and V cavities
- Fissure sealant
- Sealing hypersensitive areas
- Repair of (indirect) aesthetic restorations, temporary crown & bridge, defect margins when margins are in enamel
- Blocking out undercut
- Liner or base
- Realisation of crowns & bridges, inlays and veneers using the indirect technique in combination with GRADIA or GRADIA PLUS components (please refer to their dedicated IFU).
- Splinting of teeth in combination with fibres such as GC everStick fibres. Follow manufacturer's instructions.

CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known methacrylate monomer and methacrylate polymer allergies.

OPERATION

1. Preparations

- Hold the syringe upright and remove the wing cap by turning it anticlockwise.
- Promptly and securely attach the dispensing tip to the syringe by turning clockwise.
- Place the light protective cover until use.

Note:

Take care not to attach the dispensing tip too tightly. This may damage its screw. In order to ensure a tight connection make sure that the threads are free of remnant material.


Fig.1

2. Shade Selection

Select shade from 16 shades of XBW, BW, A1, A2, A3, A3.5, A4, B1, B2, CV, CVD, AO1, AO2, AO3, JE, AE, A, B, C, D shades are based on Vita® Shades.

*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

3. Cavity Preparation

Prepare cavity using standard techniques. Dry by gently blowing with oil free air. Rubber dam is recommended to isolate the prepared tooth from contamination with saliva, blood or sulcus fluid.

Note:

For pulp capping, use calcium hydroxide.

4. Bonding Treatment

For bonding G-aenial Universal Injectable to enamel and/or dentin, use a light-cured bonding system such as GC G-Premio BOND, GC G-aenial Bond or GC G-BOND. Follow manufacturer's instructions.

5. Placement of G-aenial Universal Injectable

- Remove the cover from the dispensing tip.
Note:
Prior to syringe material into the cavity, make sure to check the following:
1. Dispensing tip is securely attached to the syringe.
2. Any possible air from the dispensing tip is removed by gently pushing forward the syringe plunger with the tip pointing upwards until material reaches the opening of the tip (Fig. 2).


Fig.2

- Place the dispensing tip as close as possible to the prepared cavity, and slowly push the plunger to extrude material. Alternatively, displace the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.
- Place the dispensing tip as close as possible to the prepared cavity, and slowly push the plunger to extrude material. Alternatively, displace the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

Clinical Hints

In order to inject effectively, use the surface tension of the material to ensure uniformity across the entire surface of the restoration during build up. Once the required amount has been injected, release the pressure on the plunger and withdraw the syringe in a direction perpendicular to the surface. This will allow the material to separate from the dispensing tip and provides a smooth surface over the restoration.

6. Light-curing

Light-cure G-aenial Universal Injectable using a light curing unit. Keep light guide as close as possible to the surface. Refer to the following chart for Irradiation Time and Effective Depth of Cure.

Irradiation Time and Effective Depth of Cure (ISO 4049)	
10 sec. (High power LED) (> 1200mW/cm²)	20 sec. (Halogen/LED) (> 700mW/cm²)
A1, A2, A3, B1, B2, JE, AE	2.5mm
XBW, BW, A3.5, A4, AO1, AO2, CV, CVD	2.0mm
AO3	1.5mm

Note:

- Material should be placed and light-cured in layers. For maximum layer thickness, please consult above table.
- Lower light intensity may cause insufficient curing or discoloration of the material.
- For the indirect technique light-cure each layer for 30 seconds with GC LABOLIGHT LV-III, II or for 10 seconds with GC STEPLIGHT SL-I. Final cure for 3 minutes using GC LABOLIGHT LV-III.

7. Finishing and Polishing

Finish and polish using standard techniques.

SHADES

XBW, BW, A1, A2, A3, A3.5, A4, B1, B2, CV, CVD, AO1, AO2, AO3, JE, AE

STORAGE

Recommended for optimal performance, store in a cool and dark place (4 - 25°C) away from high temperatures or direct sunlight. (Shelf life: 3 years from date of manufacture)

PACKAGES

- Syringe 1.7g (1.0mL) x1, Dispensing Tip Long Needle x10, Light-protective cover x1
- Dispensing Tip Refill: Dispensing Tip Long Needle x30, Light-protective cover x3

CAUTION

- In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.

GC

CE0086

MANUFACTURED BY
GC DENTAL PRODUCTS CORP.
2-285 Torimatsu-cho, Kasugai, Aichi 486-0844, Japan

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240,
Interleuvenlaan 33, B-3001 Leuven, Belgium
TEL: +32 16 74 10 00

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