



# microSystem® - Dispenser MKII

The microSystem Dispenser with 25ml cartridges offers significant benefits over the 50ml systems currently available.

**31 %** less waste – saves money

**28 %** less weight – gentle and precise application

**26 %** smaller profile – does not restrict the intraoral view

# No blocked cartridges.


## No cross contamination.

The patented Coltène microSystem cartridges and mixing tips avoid contamination. Outflow apertures are set apart and are combined with a safety reservoir which avoids undesirable cross-contamination of material flowing back. A premature contact of both materials and subsequent undesirable polymerisation is prevented with the mixing tips.

## 31% less waste – saves money


microSystem minimises waste material with each impression with its specially developed mixing tips.

### Conventional cartridges



Conventional cartridges

### microSystem® cartridges




microSystem™ cartridges

# Gentle and precise application.

- The ergonomic microSystem Dispenser is perfectly suited to both small and large hands.
- The accurate and smooth extrusion mechanism facilitates precise intra-oral application of the impression material
- The Dispenser is easy to clean and autoclavable for complete infection control.

The picture below clearly shows the difference in size; the real tactile advantage will be immediately apparent during use.



## The perfect combination with AFFINIS®

- The compact size of the Coltène® microSystem Dispenser provides relaxed, comfortable, impression taking.
- Using microSystem you can choose AFFINIS light and regular body. The “wash” impression materials are compatible with AFFINIS heavy body, MonoBody and putties.
- The unique surface affinity and the optimized characteristics of AFFINIS light- and regular body, in combination with the microSystem Dispenser MKII, are perfectly adjusted to each other.



**AFFINIS® influences your success. All down the line.**



## Order Information

<b>6500</b>	<b>AFFINIS light body microSystem</b> 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tips (small), 20 pcs Oral Tips
<b>6600</b>	<b>AFFINIS fast light body microSystem</b> 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tips (small), 20 pcs Oral Tips
<b>6510</b>	<b>AFFINIS regular body microSystem</b> 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tips (medium), 20 pcs Oral Tips
<b>6610</b>	<b>AFFINIS fast regular body microSystem</b> 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tips (medium), 20 pcs Oral Tips

<b>6779</b>	<b>AFFINIS®PRECIOUS light body microSystem</b> 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tip Green, 20 pcs Oral Tips
<b>6780</b>	<b>AFFINIS®PRECIOUS regular body microSystem</b> 4 × 25 ml Cartridges microSystem, 20 pcs Mixing Tip Clear, 20 pcs Oral Tips

### Accessories: microSystem®

<b>6270</b>	Dispenser microSystem	1 pcs
<b>6214</b>	Mixing Tip Clear	40 pcs
<b>6215</b>	Mixing Tip Clear	100 pcs
<b>6202</b>	Mixing Tip Green	40 pcs
<b>6205</b>	Mixing Tip Green	100 pcs
<b>6210</b>	Oral Tips	100 pcs

<b>6530</b>	<b>AFFINIS putty soft</b> Base and Catalyst, each 300 ml
<b>6531</b>	<b>AFFINIS fast putty soft</b> Base and Catalyst, each 300 ml
<b>6535</b>	<b>AFFINIS putty super soft</b> Base and Catalyst, each 300 ml
<b>6520</b>	<b>AFFINIS heavy body 75 ml</b> 2 × 75 ml Cartridges, 8 pcs Mixing Tips (turquoise/large)
<b>6620</b>	<b>AFFINIS fast heavy body 75 ml</b> 2 × 75 ml Cartridges, 8 pcs Mixing Tips (turquoise/large)
<b>6547</b>	<b>AFFINIS MonoBody 75 ml</b> 2 × 75 ml Cartridges, 8 pcs Mixing Tips (turquoise/large)

### Accessories: 50/75 System

<b>6900</b>	Dispenser 50/75 ml	1 pc
<b>6552</b>	Mixing Tip Turquoise	40 pcs

<b>6487</b>	<b>AFFINIS System 360 heavy body</b> 2 × 300 ml Refill package Base, 2 × 62 ml Refill package Catalyst
<b>6485</b>	<b>AFFINIS System 360 heavy body Starter Kit</b> 1 × 300 ml Base, 1 × 62 ml Catalyst, 1 × Supporting Cartridges, 10 × Dynamic Mixing Tips
<b>6482</b>	<b>AFFINIS MonoBody System</b> 2 × 300 ml Refill package Base, 2 × 62 ml Refill package Catalyst
<b>6480</b>	<b>AFFINIS System 360 MonoBody Starter Kit</b> 1 × 300 ml Base, 1 × 62 ml Catalyst, 1 × Supporting Cartridges, 10 × Dynamic Mixing Tips

### Accessories: System 360

<b>6160</b>	Supporting Cartridges	1 pc
<b>6161</b>	Dynamic Mixing Tip Green	40 pcs

www.coltenewhaledent.com  
© 2010 Coltène/Whaledent

**Coltène/Whaledent Ltd.**  
The President Suite-A  
Kendal House, Burgess Hill  
Victoria Way  
West Sussex, RH15 9NF/U.K.  
Tel. +44 (0)1444 235486  
Fax +44 (0)1444 870640  
info@coltenewhaledent.co.uk

**Coltène/Whaledent AG**  
Feldwiesenstrasse 20  
9450 Altstätten/Switzerland  
Tel. +41 (0)71 757 53 00  
Fax +41 (0)71 757 53 01  
info@coltenewhaledent.ch

**Coltène/Whaledent Inc.**  
235 Ascot Parkway  
Cuyahoga Falls, Ohio 44223/USA  
Tel. +1 330 916 8800  
Fax +1 330 916 7077  
info@coltenewhaledent.com

**coltène  
whaledent**