

Premium-Plus-Class

The autoclaves of our new **Premium-Plus-Class** are still the **benchmark for all class B autoclaves**.

With their **twin-chamber-technology**, the biggest **colour-touch-display** and the possibility of releasing sterilization and evaluating the batch-control-system on the autoclave's screen, they set new standards on the market.

A class B cycle takes 19 minutes to sterilize and dry wrapped instruments and **only 11 minutes for unwrapped instruments**.

The clearance, labelling and traceability of wrapped sterilized equipment is facilitated with the barcode label printer MELAprint 60 which can be connected directly to autoclaves of the Premium-Plus-Class (for further information, see page 28: Documentation).

**9 kg
maximum
load!**

Device design: **tap-water connection**.

Connection is as easy as with a washing machine.

Vacuklav®40 B+

Chamber: **depth 35 cm, ø 25 cm, 18 litres.**

Tray (device holds 6 trays, dim. 29 x 19cm)

Order no.

Euro

10060

00280

Vacuklav®44 B+

Chamber: **depth 45 cm, ø 25 cm, 22 litres.**

Tray (device holds 6 trays, dim. 42 x 19cm)

Order no.

Euro

10064

00230

Device design: **stand alone**.

The only requirement for operation is an electrical socket.

Vacuklav®41 B+

Chamber: **depth 35 cm, ø 25 cm, 18 litres.**

Tray (device holds 6 trays, dim. 29 x 19cm)

Order no.

Euro

10061

00280

Vacuklav®43 B+

Chamber: **depth 45 cm, ø 25 cm, 22 litres.**

Tray (device holds 6 trays, dim. 42 x 19cm)

Order no.

Euro

10063

00230



Flex-Display As an option, the Premium-Class can also be equipped with the Flex-Display instead of the standard display that is fixed with the housing of the autoclave. In this case the autoclave is connected via a 1.5 metre long cable to the display, which can be mounted on the wall or be set up on or near by the autoclave. The Flex-Display allows you to install the autoclave in a space only 50cm high because the wall cabinet.

Please note on your order for the autoclave: "With Flex-Display".

Euro

Upgrading the Premium-Class autoclave from the black/white display **to the new colour-touch-display** is possible with little effort. The scope of delivery includes the colour-touch display, clearance software and the components required for the conversion.

The price does not include the costs of the approx. 1 hour of working time and any travel expenses.

Order no.

Euro

67205

Three different mounts for your choice

Premium-Class autoclaves are fitted with the mount >C< for holding trays, MELAstore-Boxes or sterilization containers as standard. Should you require another mount than the standard, please inform us with the order if you like mount >D< (for high cassettes, e.g. MELAstore-Box 200) or mount >E< for 6 standard trays and 2 additional narrow trays (to be ordered separately) for an optimum use of the chamber.



Mount >C<
for 6 trays



Mount >C<
(turned 90°) for
3 MELAstore-Box 100



Mount >D<
for 2 high cassettes,
e.g. MELAstore-Box 200



Mount >E<
for 6 standard trays
and 2 narrow trays

	Order no.	Euro
For chamber depth 35cm:		
Mount >C<	40232	
Mount >D<	46840	
Mount >E<	40237	
Narrow tray (short) for Mount >E<	01320	

	Order no.	Euro
For chamber depth 45cm:		
Mount >C<	40242	
Mount >D<	46840	
Mount >E<	40238	
Narrow tray (long) for Mount >E<	01310	

MELAstore® Trays • Wash trays for the automatic reprocessing

The reprocessing of instruments with the **integrated system** of MELAstore-Trays and MELAstore-Boxes **saves time and money**. With MELAstore you **simplify your workflow** of the instruments from the patient treatment to the reprocessing and finally the storage of sterilized goods.

The MELAstore-Trays are made of **high-quality stainless steel** and are optimized for the loading of MELAtherm. The wide opened mesh of the MELAstore-Trays ensure **best cleaning results**. For the comprehensive and complete documentation of the instrument reprocessing, the MELAstore-Trays can be equipped with identification plates.

MELAstore-Trays offer:

- significant time saving in the loading and unloading of the instruments from MELAtherm
- increasing the quality of the overall reprocessing through standardization and documentation
- a safe handling of the instruments by fixing them into silicone bars
- protection against personal injury of the practice team



MELAstore Tray 50



Silicone bar of MELAstore-Tray 50



MELAstore Tray 100



Silicone bar of MELAstore-Tray 100



MELAstore Tray 200



Silicone bar of MELAstore-Tray 200

	Order-no.	Euro
MELAstore®-Tray 50 Delivered incl. two silicone bars for up to 8 instruments with a length of max. 17.5cm. Width: 11.8cm x Length: 18cm x Height: 3cm	01180	
MELAstore Silicone bar of MELAstore-Tray 50	82960	
MELAstore®-Tray 100 Delivered incl. three silicone bars for up to 14 instruments with a length of max. 27cm. Width: 17.6cm x Length: 27.5cm x Height: 3cm	01181	
MELAstore Silicone bar of MELAstore-Tray 100	82970	
MELAstore®-Tray 200 Delivered incl. three silicone bars for up to 13 instruments with a length of max. 27cm. Width: 17.6cm x Length: 27.5cm x Height: 4.3cm	01182	
MELAstore Silicone bar of MELAstore-Tray 200	82980	

MELAstore® Box • sterilization containers according DIN 58953-9

For sterilization, the MELAstore-Trays are inserted in the corresponding MELAstore-Box. The **lid** of the MELAstore-Boxes are **perforated** and equipped with a **holder for a filter**. The filters are available as **disposable paper filters** or as **permanent filters** and have to be ordered separately. For marking the MELAstore-Boxes MELAG offers identification plates with or without barcode that are fixed in the holder of the MELAstore-Box. There are also disposable seals available, with or without process indicator.



MELAstore Box 100



MELAstore Box 200

	Order-no.	Euro
MELAstore®-Box 100 For the acceptance of one MELAstore-Tray 100 or of 2 MELAstore-Trays 50 Width: 19cm x Length: 31.2cm x Height: 4.6cm	01191	
MELAstore permanent filter , 2 pcs. (for up to 1,000 cycles)	01193	
MELAstore disposable paper filter , 100 pcs.	01194	

	Order-no.	Euro
MELAstore®-Box 200 For the acceptance of one MELAstore-Tray 200 or of 2 MELAstore-Trays 100 Width: 19cm x Length: 31.2cm x Height: 6.5cm	01192	
MELAstore permanent filter , 2 pcs. (for up to 1,000 cycles)	01193	
MELAstore disposable paper filter , 100 pcs.	01194	

MELAstore® permanent filters • according European standard EN 868-2

MELAstore permanent filters are designed for the use in the MELAstore-Box 100 and in the MELAstore-Box 200. In case of damages or soiling the filter has to be replaced, latest after 1,000 sterilization cycles.

MELAstore® permanent filter

for use with MELAstore-Box 100 and MELAstore-Box 200
(Packaging unit: 2 pieces)

Order no.	Euro
01193	



MELAstore permanent filter

MELAstore® disposable paper filters • according European standard EN 868-2

The MELAstore disposable paper filters are designed for the use in the MELAstore-Box 100 and in the MELAstore-Box 200. For each sterilization cycle a new paper filter is inserted into the lid of the MELAstore-Box.

MELAstore® disposable paper filter

for use with MELAstore-Box 100 and MELAstore-Box 200
(Packaging unit: 100 pieces)

Order no.	Euro
01194	



MELAstore disposable paper filter

MELAstore® disposable seal • with or without process indicator

For an easy and safe sealing of the MELAstore-Boxes and the MELAG sterilization containers. The disposable seals are available with or without a process indicator.

MELAstore® disposable seal

without process indicator. Packaging unit: 100 pieces

Order no.	Euro
01198	

MELAstore® disposable seal

with process indicator. Packaging unit: 100 pieces

01199	
-------	--

Disposable seal without indicator



Disposable seal with indicator



MELAstore® holder for identification plate for MELAstore-Tray

This holder is fixed at the lid of a MELAstore-Tray. The identification plate (order no. 01195) can be inserted into this holder.

MELAstore® holder for identification plate

for fitting on the MELAstore-Tray 50, 100 and 200

Order no.	Euro
01197	



Holder for identification plate

MELAstore® identification plate MELAstore-Tray / MELAstore-Box

This identification plate can be printed with a text of up to 15 characters and / or additionally with a barcode. The identification plates will be printed as requested by the practice.

The identification plate will be inserted into the holder of a MELAstore-Box or into the holder for an identification plate (order no. 01197) at the MELAstore-Tray.

MELAstore® identification plate

incl. printing of up to 15 characters and additionally barcode

Order no.	Euro
01195	



Identification plate with characters and barcode

MELAstore® identification plate tall for MELAstore®-Box

The extra-tall identification plate (not labelled) is for affixing the labels of a label-printer, e. g. the MELAprint 60. This identification plate will be inserted into the holder of a MELAstore-Box.

MELAstore® identification plate tall

for use with MELAstore-Box 100 and MELAstore-Box 200
(Packaging unit: 10 pieces)

Order no.	Euro
01196	



Identification plate tall with label of MELAprint 60

MELAfol® • See-through sterilization package according EN 868-5

Available in rolls, or as individual bags. With an indicator that turns from rose to brown during sterilization.

Reels (Width x Length)	Order no.	Euro
MELAfol®502 (5cm x 200m)	00502	
MELAfol®752 (7.5cm x 200m)	00752	
MELAfol®1002 (10cm x 200m)	01002	
MELAfol®1502 (15cm x 200m)	01502	
MELAfol®2002 (20cm x 200m)	02002	
MELAfol®2502 (25cm x 200m)	02502	

Pouches (Width x Length)	Order no.	Euro
MELAfol®501 (5cm x 25cm), 1000 Stück	00501	
MELAfol®751 (7.5cm x 25cm), 1000 Stück	00751	
MELAfol®1001 (10cm x 25cm), 1000 Stück	01001	
Pouches with fold (Width x Length)	02051	
MELAfol®2051 (20cm x 50cm), 100 pieces for sealing standard-tray-cassettes		

Foil holder

For sterilizing instruments sealed in MELAfol upright for better drying results.

	Order no.	Euro
For autoclaves 29 VS+ / 30 B+ / 31 B+ / 40 B+ / 41 B+	22410	
For autoclaves 23 VS+ / 23 B+ / 24 B+ / 43 B+ / 44 B+	22420	
For autoclaves MELAtronic 23EN	00283	
For autoclaves MELAtronic 15EN+	70670	



MELAfol reels and bags



Foil holder

MELAG seal check • routine control of sealing devices

For the weekly verification of the quality of the seal seam.

	Order no.	Euro
MELAG seal check , 100 test strips	01079	



MELAG seal check

Sterilization containers according to European standard EN 868-8

The MELAG sterilization containers consist of a lower part made of anodized aluminium and an upper part of stainless steel. The base and lid of the sterilization containers are perforated and fitted with single-use paper filters. 100 single-use paper filters are included in the scope of delivery.

Type	Outer dimensions (Length x Width x Height, cm)	Order no.	Euro
15K	18 x 12 x 4.5	01151	
15M	35 x 12 x 4.5	01152	
15G	35 x 12 x 8.0	01153	
17K	20 x 14 x 5.0	01171	
17M	41 x 14 x 5.0	01172	
17G	41 x 14 x 9.0	01173	
23M	42 x 16 x 6.0	01231	
23G	42 x 16 x 12.0	01232	
28M	32 x 16 x 6.0	01284	
28G	32 x 16 x 12.0	01285	



Type K

Type M

Type G

Paper filter for sterilization containers according to European standard EN 868-2

A new paper filter for the base and lid of the sterilization container are required for each sterilization cycle.

Single-use paper filters (100 pieces)	Dimensions (Length x Width, cm)	Order no.	Euro
For sterilization container 15K	11.5 x 10.5	27325	
For sterilization container 15M / 15G	29.5 x 10.5	27335	
For sterilization container 17 K	13.0 x 12.0	27345	
For sterilization container 17M / 17G	35.5 x 12.5	27355	
For sterilization container 23M / 23G	35.2 x 13.2	26515	
For sterilization container 28M / 28G	24.2 x 14.0	27365	



Single-use paper filters, 100 pieces