

W Add "W" to part number for individually packaged arch wires where appropriate.

MAXILLARY	Form I		Form III**	
		Dimpled		Dimpled
	.012	857-502 W	857-573 W	857-702I W
	.014	857-504 W	857-574 W	857-704I W
	.016	857-506 W	857-576 W	857-706I W
	.018	857-508 W	857-578 W	857-708I W
	.020	857-510 W	857-579 W	857-710I W
	.016 x .016	857-540 W	857-540D W	857-731I W
	.016 x .022	857-541 W	857-541D W	857-732I W
	.017 x .025	857-537 W	857-537D W	857-733I W
	.018 x .018	857-545 W	857-545D W	857-734I W
	.018 x .025	857-542 W	857-542D W	857-735I W
	.019 x .025	857-543 W	857-543D W	857-736I W

MANDIBULAR	Form I		Form III**	
		Dimpled		Dimpled
	.012	857-512 W	857-583 W	857-712I W
	.014	857-514 W	857-584 W	857-714I W
	.016	857-516 W	857-586 W	857-716I W
	.018	857-518 W	857-588 W	857-718I W
	.020	857-520 W	857-589 W	857-720I W
	.016 x .016	857-550 W	857-550D W	857-741I W
	.016 x .022	857-551 W	857-551D W	857-742I W
	.017 x .025	857-547 W	857-547D W	857-743I W
	.018 x .018	857-555 W	857-555D W	857-744I W
	.018 x .025	857-552 W	857-552D W	857-745I W
	.019 x .025	857-553 W	857-553D W	857-746I W



3524 Washington Avenue
Sheboygan, Wisconsin, 53081 USA

Phone: 1 920 457 5051
Toll Free: 1 800 558 7687
Email: info@americanortho.com
Web: americanortho.com




©2016 American Orthodontics Corporation
All Rights Reserved

FOLLOW US ON
f y in t

NT3™

SE NiTi Arch Wire



NT3
SE NiTi Arch Wire

SUPERELASTIC
SHAPE MEMORY
SPRINGBACK

American Orthodontics' NT3 Nickel Titanium Arch Wire is the latest in AO's extensive arch wire line. Manufactured at AO's state-of-the-art manufacturing facility in Sheboygan, Wisconsin, this high-quality wire delivers exceptional superelasticity, shape memory, and springback for precise control and predictable results.

Available arch form options include VLP* Arch Form, Natural Arch Form I, and Natural Arch Form III**.

- Manufactured using premium Japanese NiTi material
- Highly resistant to permanent deformation
- Polished surface finish
- Consistent loading and unloading forces
- Thoroughly tested to ensure consistent forces between lots
- Sterilizable packaging options
- Pre-loaded stops and dimple options
- 10 per package

* Compare VLP to the Damon® Arch Form, a registered trademark of Ormco
** Compare Form III to the European Arch Form

W

Add "W" to part number for individually packaged arch wires where appropriate.

UNIVERSAL	VLP*		
		Dimpled	Pre-loaded stops
	.012	857-112 W	857-132 W 857-112PW
	.014	857-114 W	857-134 W 857-114PW
	.016	857-116 W	857-136 W 857-116PW
	.018	857-118 W	857-138 W 857-118PW
	.014 x .025	857-120 W	857-140 W 857-120PW
	.016 x .025	857-121 W	857-141 W 857-121PW
	.018 x .018	857-122 W	857-142 W 857-122PW
	.018 x .025	857-123 W	857-143 W 857-123PW
	.019 x .025	857-124 W	857-144 W 857-124PW
	.020 x .020	857-125 W	857-145 W 857-125PW

Lower VLP* Options	
	Small Dimple
.012	857-132LSP W
.014	857-134LSP W
.016	857-136LSP W
.018	857-138LSP W
.014 x .025	857-140LSP W
.016 x .025	857-141LSP W
.018 x .018	857-142LSP W
.018 x .025	857-143LSP W
.019 x .025	857-144LSP W
.020 x .020	857-145LSP W

AVAILABLE ARCH FORMS

NATURAL ARCH FORM I

Well suited for pre-adjusted appliances with programmed In/Out relationships, Form I can be used throughout treatment with little or no manipulation. Its moderately expanded anterior segment provides a slight over correction to allow for settling after removal.

NATURAL ARCH FORM III** - European Arch Form

Form III is slightly narrower than Form I in the anterior region, yet is well suited for preadjusted appliances. Due to the gentle curve of the arch form in the anterior segment, Form III is ideal for single-wing bracket systems as well as self ligating bracket systems where low friction is critical to appliance performance.

VLP* ARCH FORM

The VLP arch form was designed to maximize the benefits of a passive self ligating bracket system. The wider arch width in the posterior allows for full expression of the expansion characteristics of low friction passive systems.

