level indicator O/KULZER 3 000 Dynamix*speed 2 -

Dynamix[®]*speed 2* Precision at the push of a button – now even better

- **>>** For your convenience: Additional markings on the level indicator allow you to assess whether the remaining quantity is sufficient for another impression (upper/lower jaw). This helps to improve the impression-taking process in the practice.
- **>>** For the environment: Avoids waste as you can make the best possible use of the contents of the cartridge.

KULZER

Flexitime[®] Product Sheet

Giving a hand to oral health.

New material

Dynamix[®]speed 2


Precision at the push of a button

When it comes to high-speed mixing, time-saving and simultaneously obtaining a consistent, void-free mixture, the decision is easily made: practical cartridges and the functional Dynamix®speed 2 mixer are perfectly matched in this well-thought-out system featuring all the advantages of automatic mixing for impression materials.

What's new with Dynamix®speed 2:

- For your convenience: Additional markings on the level indicator allow you to assess whether the remaining quantity is sufficient for another impression (upper/lower jaw). This helps to improve the impression-taking process in the practice.
- For the environment: Avoids waste as you can make the best possible use of the contents of the cartridge.

New material



66093755	DYNAMIX SPEED 2 EU
40004353	Dynamix Mixing Tips (50 pcs)
66047385	Dynamix Fixation Ring Yellow (2 pcs)
66044882	Flexitime Fast & Scan Dynamix Putty Refill
66045755	Flexitime Fast & Scan Dynamix Heavy Tray Refill
40003992	Flexitime Dynamix Putty Refill
40003991	Flexitime Dynamix Heavy Tray Refill
40003990	Flexitime Dynamix Monophase Refill

Reliable in quality and precision

 Reproducible material properties. Consistent quality of mix, independent of operator, for a constant precision of fit.

including for Putty viscosities. Easy and reliable handling. Savings in time and material. 3 selectable velocities.

restorations.

High-speed efficiency

Art.-Code