

# LMActivator™

LM

feel the  
difference



## A million smiles found, and more to come.

### Meet the new **My LMActivator™**

A renowned and clinically proven appliance with two decades of history in early orthodontic treatment, is now available as a new and enhanced product line: meet My LM-Activator™ by LM-Dental™.

My LM-Activator™ is a **prefabricated silicone aligner** for growing happy smiles: easy to use, gently supporting healthy and natural growth of the face and jaw. The ideal time for treatment is in early mixed dentition, as research indicates that malocclusions in deciduous dentition tend to become more severe as the dentition develops. By intervening and taking action to correct malocclusions early, the need for further treatment can be significantly reduced and even eliminated.

Developed in close collaboration with clinicians, the new My LM-Activator™ product line has been **enhanced with optimized dimensions and advanced features** for user-friendliness and comfort.

My LM-Activator™ represents a leap forward in our dedication to providing high quality orthodontic solutions for growing patients.



NEW!



Introducing the new, enhanced  
**My LM-Activator™** silicone activator and  
aligner for early orthodontic treatment.

## My LMActivator™

# What's new in it?

Discover enhanced features in My LM-Activator™




## Streamlined anatomy

My LM-Activator™ presents an overall streamlined and light-weight design **improving fit and comfort for the patient**. It allows more space for the canines and for development of the upper arch and jaw.



## Harmonized design and sizing

A harmonized design over the size range makes it even easier to find a good fit for your patient. The ultimate **high-gloss finish** is combined with the **superior wear comfort** that LM-Activator™ is well known for.




## Optimized vertical control

The **improved thickness profile of the occlusal plane** enables vertical control, empowering the development of optimal occlusal conditions in mixed dentition.

# Seamless transition for current LM-Activator™ users

My LM-Activator™ has been developed in **close cooperation with existing LM-Activator™ users and orthodontic specialists** with long clinical experience of LM-Activator™ use. All design and feature improvements are based on the needs and insights of these experts.

As the treatment principles remain unchanged, current LM-Activator™ users can seamlessly transition to the My LM-Activator™. The new product line includes all well-established and the most commonly-used models and sizes, enabling a gradual simplification of the LM-Activator™ product assortment.



## My LM-Activator™ product range

My LM-Activator™ is the new product line of LM-Activator™. With streamlined and lightweight design combined with further optimized features and morphology it is continuing the established path of LM-Activator™. Experience both enhanced comfort and optimized fit in the new My LM-Activator™. The most common models and sizes for your patient are easy to find.



My LM-Activator™ Low Short ●

| SIZE | CODE    |
|------|---------|
| 35   | 94435LS |
| 40   | 94440LS |
| 45   | 94445LS |
| 50   | 94450LS |
| 55   | 94455LS |
| 60   | 94460LS |
| 65   | 94465LS |
| 70   | 94470LS |

Model with a shorter molar section when the molars have not yet erupted.

Delivered in a yellow box.



My LM-Activator™ Low Long ●

| SIZE | CODE    |
|------|---------|
| 40   | 94440LL |
| 45   | 94445LL |
| 50   | 94450LL |
| 55   | 94455LL |
| 60   | 94460LL |
| 65   | 94465LL |

Model with a longer molar section when the second molars have erupted.

Delivered in a blue box.



My LM-Activator™ High Short ●

| SIZE | CODE    |
|------|---------|
| 40   | 94440HS |
| 45   | 94445HS |
| 50   | 94450HS |
| 55   | 94455HS |
| 60   | 94460HS |
| 65   | 94465HS |

Model with a shorter molar section. The High model is thicker in the region of second premolars and molars. It is designed for treating skeletal and dentoalveolar Open bite cases.

Delivered in an orange box.



My LM-Activator™ High Long ●

| SIZE | CODE    |
|------|---------|
| 45   | 94445HL |
| 50   | 94450HL |
| 55   | 94455HL |
| 60   | 94460HL |

Model with a longer molar section. The High model is thicker in the region of second premolars and molars. It is designed for treating skeletal and dentoalveolar Open bite cases.

Delivered in a green box.

More information [www.lm-activator.com](http://www.lm-activator.com)