POCKETGUIDE

ORTH**O**easy®





Pin .022" ø 1.7 mm

Color: Rosé

Length: 6 mm

Diameter: ø 1.7 mm

Thread length: 5.8 mm

Total length: 9.8 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Cross slot .022"

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1101A2306 OrthoEasy® Pin



Rosé

6 mm, ø 1.7 mm

Color: Pink

Length: 8 mm

Diameter: ø 1.7 mm

Thread length: 7.8 mm

Total length: 11.8 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Cross slot .022"

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1101A2308

OrthoEasy® Pin



Pink

8 mm, ø 1.7 mm

Color: Gold

Length: 10 mm

Diameter: ø 1.7 mm

Thread length: 9.8 mm

Total length: 13.8 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Cross slot .022"

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1101A2310 OrthoEasy® Pin



Gold 10 mm, ø 1.7 mm

Color: Silver

Length: 12 mm

Diameter: ø 1.7 mm

Thread length: 11.8 mm

Total length: 15.8 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Cross slot .022"

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1101A2312

OrthoEasy® Pin



Silver 12 mm, ø 1.7 mm



Pal ø 1.7 mm

Color: Pink

Length: 8 mm

Thread length: 7.8 mm

Total length: 12.6 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1201A2308

OrthoEasy® Pal



Pink

8 mm, ø 1.7 mm

Color: Gold

Length: 10 mm

Thread length: 9.8 mm

Total length: 14.6 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1201A2310OrthoEasy® Pal



Gold 10 mm, ø 1.7 mm

Color: Silver

Length: 12 mm

Thread length: 11.8 mm

Total length: 16.6 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1201A2312 OrthoEasy® Pal



Silver 12 mm, ø 1.7 mm

Color: Blue

Length: 14 mm

Thread length: 13.8 mm

Total length: 18.6 mm

External thread diameter: 1.7 mm

Fracture torque: \geq 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1201A2314 OrthoEasy® Pal



Blue 14 mm, ø 1.7 mm



Pal Long Neck ø 1.7 mm

Color: Pink

Length: 8 mm

Thread length: 7.8 mm

Total length: 14.7 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1221A1708

OrthoEasy® Pal Long Neck



Pink

8 mm, ø 1.7 mm

Color: Gold

Length: 10 mm

Thread length: 9.8 mm

Total length: 16.7 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1221A1710

OrthoEasy® Pal Long Neck



Gold

10 mm, ø 1.7 mm

Color: Silver

Length: 12 mm

Thread length: 11.8 mm

Total length: 18.7 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1221A1712

OrthoEasy® Pal Long Neck



Silver 12 mm, ø 1.7 mm

Color: Blue

Length: 14 mm

Thread length: 13.8 mm

Total length: 20.7 mm

External thread diameter: 1.7 mm

Fracture torque: > 30 Ncm

Insertion torque: < 20 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.0 to max. 1.3 mm

1221A1714

OrthoEasy® Pal Long Neck



Blue 14 mm, ø 1.7 mm



Pal Long Neck ø 2.3 mm

Color: Pink

Length: 8 mm

Thread length: 7.8 mm

Total length: 14.6 mm

External thread diameter: 2.3 mm

Fracture torque: > 70 Ncm

Insertion torque: < 40 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.7 mm

1221A2308

OrthoEasy® Pal Long Neck



Pink 8 mm, ø 2.3 mm

Color: Gold

Length: 10 mm

Thread length: 9.8 mm

Total length: 16.6 mm

External thread diameter: 2.3 mm

Fracture torque: > 70 Ncm

Insertion torque: < 40 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.7 mm

1221A2310

OrthoEasy® Pal Long Neck



Gold 10 mm, ø 2.3 mm

Color: Silver

Length: 12 mm

Thread length: 11.8 mm

Total length: 18.6 mm

External thread diameter: 2.3 mm

Fracture torque: > 70 Ncm

Insertion torque: < 40 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.7 mm

1221A2312

OrthoEasy® Pal Long Neck



Silver 12 mm, ø 2.3 mm

Color: Blue

Length: 14 mm

Thread length: 13.8 mm

Total length: 20.6 mm

External thread diameter: 2.3 mm

Fracture torque: > 70 Ncm

Insertion torque: < 40 Ncm

Head geometry: Internal thread M1.7

Pilot hole – ø recommended: 1.7 mm

1221A2314

OrthoEasy® Pal Long Neck



Blue 14 mm, ø 2.3 mm

Thread length: 6 mm

Variants: 1101A2306 OrthoEasy® Pin



1101A2306

Thread length: 8 mm

Variants:

1101A2308 OrthoEasy® Pin OrthoEasy® Pal 1201A2308

1221A1708

OrthoEasy® Pal Long Neck 1.7 mm OrthoEasy® Pal Long Neck 2.3 mm 1221A2308

1101A2308

1201A2308





1221A1708



1221A2308

Thread length: 10 mm

Variants:

1101A2310 OrthoEasy® Pin
1201A2310 OrthoEasy® Pal
1221A1710 OrthoEasy® Pal Long Neck 1.7 mm
1221A2310 OrthoEasy® Pal Long Neck 2.3 mm



1101A2310

- Section .

1201A2310



1221A1710



1221A2310

Thread length: 12 mm

Variants:

1101A2312 OrthoEasy® Pin1201A2312 OrthoEasy® Pal

1221A1712 OrthoEasy® Pal Long Neck 1.7 mm

1221A2312 OrthoEasy® Pal Long Neck 2.3 mm

1

1101A2312



1201A2312



1221A1712



1221A2312

Thread length: 14 mm

Variants: 1201A2314 1221A1714 1221A2314 1201A2314 OrthoEasy® Pal OrthoEasy® Pal Long Neck 1.7 mm 1221A1714 OrthoEasy® Pal Long Neck 2.3 mm 1221A2314

Abutments for all Pal

High profile with rotation lock
With hook and slot
Suitable for two elastics

Slot: .022" x .028"

Space for elastics below hook: Ø 1.0 mm

Passage for elastics: 0.5 mm

Retention force: > 10 N

Height: 3.9 mm

Outer ø: 5.0 mm

1299-1015 with hook + slot



High profile without rotation lock
With tube
For easy adaption of appliances such as e.g. RPE, Frog II

Tube ø: 1.5 mm

Height: 3.9 mm

Outer ø: 5.0 mm

1299-1016 with round tube





Flat profile without rotation lock For welding

Height: 2.1 mm

Outer ø: 5.0 mm

Height of welding surface: 1.1 mm

1299-1004 flat



High profile with rotation lock For welding

Height: 3.9 mm

Outer ø: 5.0 mm

Height of welding surface: 1.9 mm

1299-1014 high





High profile, without rotation lock For welding

Height: 3.9 mm

Outer ø: 5.0 mm

Height of welding surface: 1.9 mm

1299-1018 high





High profile without rotation lock

Easy mount — with cut-out for lateral fitting of abutment

For welding

Height: 3.9 mm
Outer ø: 5.0 mm
Height of welding surface: 1.9 mm

1299-1019 high, with cut-out





Fixation plates for all Pal

Fixation plate with parallel wire

For appliances for mesialization or distalization of molars in the upper jaw

Diameter: Ø 1.1 mm
Length fixation element: 13.5 mm
Total length: 100 mm

1299-0100 parallel wire

Fixation plate with orthogonal wire

For appliances for mesialization or distalization of molars in the upper jaw

Diameter: ø 1.1 mm
Length fixation element: 13.5 mm
Total length: 100 mm

1299-0101 orthogonal wire



Fixation plate with wire in T-shape


Total length – heigth:

For appliances for mesialization or distalization of molars in the upper jaw

65 mm

Diameter: Ø 1.1 mm Length fixation element: 13.5 mm Total length – width: 100 mm

1299-0102 wire in T-shape



Insertion instruments

Blades for contra-angled handpiece:

1199-0007 1199-0008

1300-0002-010

OrthoEasy® Octagonal blade, 5 mm

OrthoEasy® Octagonal blade, 10 mm

OrthoEasy® Blade for insertions guides

(with depht stop)



1199-0007



1199-0008



1300-0002-010

Blades for manual driver:

1199-0006 OrthoEasy® Octagonal blade, 10 mm 1199-00015 OrthoEasy® Octagonal blade, 45 mm



1199-0002: OrthoEasy® handle for screw driver



1199-0006



1199-0015

Accessories for all Pal

Retaining screw for fixation of abutments or fixation plates



1299-1006: Blade for retaining screw with saftey hold



1299-0007: SmartDrive torque driver
with torque limiting for blade
for retaining screw and
SmartJet™ blade

1299-1003Retaining screw





Distal-mesial slider

For appliances for mesialization or distalization of molars in the upper jaw



322-0110SmartJetTM Tube



Accessories for Lab

For the preparation of the appliances in the laboratory



1299-1006: Blade for retaining screw with saftey hold



1299-0007: SmartDrive torque driver
with torque limiting for blade
for retaining screw and
SmartJet™ blade

1299-1001Lab analog



FORESTADENT

Bernhard Förster GmbH Westliche Karl-Friedrich-Str. 151 75172 Pforzheim Phone: +49 (0) 7231 459-0 Fax: +49 (0) 7231 459-102 www.forestadent.com

info@forestadent.com

