



OptiBond™ Universal 360

Single Component Universal Adhesive

All Techniques. Any Material. One Bonding Solution. The One for All.

Meet OptiBond Universal 360, the revolutionary new bonding agent that offers the best outcome to every clinician, regardless of what materials they use or which techniques they follow. No matter what your goals are, OptiBond Universal 360 can help your dental practice meet them with ease.



Universal, Multi-Purpose Solution

Powered by Kerr's patented GPDM technology, now enriched with MDP, OptiBond Universal 360 is a truly universal system offering reliable bond strength to all surfaces and substrates, with any resin cement, even in the absence of light. This one-bottle solution offers the flexibility to be used with any etching technique in both direct and indirect restorations.



Exceptional Bond Strength

OptiBond Universal 360 offers equal or superior bond strength results in combination with Kerr and leading competitor cements, even in total dark cure mode, without auxiliary products.



Improved Bottle Delivery

Thanks to the all new ergonomic bottle design, you'll only need one hand for opening, dispensing, and closing. Plus, you can use the dropper to deliver a precise drop size at every application, without worrying about leakage before and after. This leads to less waste and more efficiency.

Revolutionize Your Practice with One Bottle



Improved dark cure capabilities

OptiBond Universal 360 optimizes the balance between hydrophilicity and hydrophobicity in its formula to prevent phase separation at the adhesive and resin cement interface. This results in improved complete polymerization, even in total dark cure mode.



Three Solvents in One

Kerr's patented Ternary Solvent System technology and optimal hydrophobicity and hydrophilicity balance allow OptiBond Universal 360 to deliver excellent penetration into dentinal tubules, giving you extraordinary bond strength and protection against microleakage and post-op sensitivity.



Easy Storage


Store OptiBond Universal 360 at an ambient temperature. Unlike other bonding agents, it doesn't require refrigeration.



Minimal Post-Op Sensitivity

OptiBond Universal 360 ensures optimal tooth penetration to help you reduce post-op sensitivity risk and increase patient satisfaction.

THE BENEFIT OF COMBINING TWO MONOMERS

PROPERTY	 GPDM	 MDP	 GPDM + MDP
Acidity	Lower pH.	Higher pH.	OptiBond Universal 360 Adhesive formula with optimized hydrophilicity-hydrophobicity and pH, for wider range of substrate compatibility and convenient storage requirement.
Hydrophilicity/ Hydrophobicity	More hydrophilic.	More hydrophobic.	
Adhesion anchoring group	Adhesion via phosphate group.	Adhesion via phosphate group.	
Crosslinking to resin	Two functional groups for crosslinking to resin.	One functional group for crosslinking to resin.	

Ordering Information

Part #	Item Description
36820	OptiBond Universal 360 - Bottle Kit Content: 1x OptiBond Universal 360 bottle 5ml, 1x Gel Etchant 3g syringe, 25x disposable mixing wells, 100x disposable applicator tips, technique guide.
36821	OptiBond Universal 360 - Unidose Kit Content: 100x OptiBond Universal 360 unidose applications 0.18 ml, 1x Gel Etchant 3g syringe, 100x disposable applicators, technique guide.
36822	OptiBond Universal 360 - Bottle Refill Content: 1x OptiBond Universal 360 bottle 5ml



For more information, visit
to.kerrdental.com/optibond360

Sweorto AB | S:t Johannesgatan 2, | 211 46 Malmö | Sweden
 Kerr Dental Nordics & Baltics | kerrdental.com

©2025 Kerr Corporation. All trademarks applicable to Kerr and Kerr products are property of Kerr Corporation. All other trademarks are property of their respective owners. All Rights Reserved. MKT-25-0143

