

# CONVINCING RETRACTION



RetraXil®

Starter Pack

☞ Adstringierende Retraktionspaste

☞ Astringent retraction paste



KULZER

MITSUI CHEMICALS GROUP

## RetraXil®

### Good sulcus exposure – minimally invasive application

RetraXil® provides a cleanly opened gingival sulcus and ensures beautifully exposed, blood-free preparation margins. The innovative twisting aid supports application – the cannula stays securely in place.

- » Atraumatic application
- » Adequate tissue management due to high viscosity.
- » Multidose stable syringe. No curls, no tear-off, non-stick.
- » Works with analogue and digital impressions.

Giving a hand to oral health.



**KULZER**  
MITSUI CHEMICALS GROUP



# RetraXil® – Convincing Retraction

Efficient, gingival retraction with highest standards in handling and safety

## Safety during application

- No risk of systemic overdose
- Blue colour contrast to the gingiva and teeth offers high visibility and good control when rinsing
- Twisting aid offers secure fit of the cannula on the syringe and enables the cannula to be twisted on and off safely

## Good handling features

- Thin cannula offers precise and atraumatic application into the sulcus and interproximal areas
- Compatible with analog and digital impression techniques

## Economical and time-saving

- Fewer work steps compared to retraction cord technique


## Dentists' first impressions about RetraXil®

"Great and strong addition to the market and loved the product. Will definitely be using it again and again."













"I was extremely pleased with the consistency of the paste. Easy to wash away leaving no particles behind and dried the sulcus efficiently."

"RetraXil® would be my choice over astringent capsules as it is more user-friendly and there isn't any residue."


Based on a 12-item survey in different countries – data on file.



**RetraXil® offers atraumatic and easy access to the sulcus thanks to the thinnest cannula on the market**

RetraXil®	Competitor 1	Competitor 2	Competitor 3	Competitor 4	Competitor 5
					
Ø 0.9 mm	Ø 1.25 mm	Ø 1.25 mm	Ø 1.4 mm	Ø 1.4 mm	Ø 1.6 mm
					
Cannula	PLT	Cannula	PLT	PLT	Cannula

Art. Code	Article	Content
<b>Products</b>		
66077159	<b>RetraXil Starter Pack</b>	4 x 1 g RetraXil 2 x Twisting aid 20 x Cannulae
66077158	<b>RetraXil Value Pack</b>	20 x 1 g RetraXil 10 x Twisting aid 100 x Cannulae
<b>Refills</b>		
66087214	<b>RetraXil Twisting Aid Refill</b>	2 x 10 Twisting Aids
66087215	<b>RetraXil Cannula Refill</b>	5 x 20 Cannulae



[Kulzer.com/Retraxil](https://www.kulzer.com/Retraxil)