

LM

Adapt. Shape. Replicate.

feel the difference



NEW!

LMArte™ Replica Posterior



REPLICATE EASILY

LM-Arte™ Replica Posterior
448-449ES

Clever instrument designed for simplifying occlusal modeling of posterior teeth in direct composite restorations. The instrument makes occlusal modeling fast, simple and predictable. Occlusal anatomy of single tooth can be made easily guided by the buccal cusps using replica technique with this specific instrument.



“ I can make modeling within seconds with simple steps by using Replica Posterior instrument and technique.”

Dr. Osama Shaalan,
Gold Member of Styleitaliano™

Scan the QR-code to watch the how-to-use -video.



Adapt.

Shape.

Replicate.



 LMDental™

LM-Instruments Oy | Norrbyn rantatie 8 | FI-21600 Parainen, Finland | info@lm-dental.com | www.lm-dental.com | www.lm-dental.com/patents

 